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Vol.2

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- 1) Hannink RHJ, Kelly PM, Muddle BC (2000) Transformation toughening in zirconia-containing ceramics. *J Am Ceram Soc* 83:461–487
- 2) Kohorst P, Dittmer MP, Borchers L, Stiesch-Scholz M (2008) Influence of cyclic fatigue in water on the load-bearing capacity of dental bridges made of zirconia. *Acta Biomater* 4:1440–1447
- 3) Kohorst P, Herzog TJ, Borchers L, Stiesch-Scholz M (2007) Load-bearing capacity of all-ceramic posterior four-unit fixed partial dentures with different zirconia frameworks. *Eur J Oral Sci* 115:161–166
- 4) Tinschert J, Natt G, Hassenpflug S, Spiekermann H (2004) Status of current CAD/CAM technology in dental medicine. *Int J Comput Dent* 7:25–45
- 5) Raigrodski AJ (2004) Contemporary materials and technologies for all-ceramic fixed partial dentures: a review of the literature. *J Prosthet Dent* 92:557–562
- 6) Denry I, Kelly JR (2008) State of the art of zirconia dental applications. *Dent Mater* 24:299–307
- 7) Rosentritt M, Behr M, Thaller C, Rudolph H, Feilzer A (2009) Fracture performance of computer-aided manufactured zirconia and alloy crowns. *Quintessence Int* 40:655–662
- 8) Comlekoglu M, Dundar M, Ozcan M, Gungor MA, Gokce B, Artunc C (2009) Influence of cervical finish line type on the marginal adaptation of zirconia ceramic crowns. *Oper Dent* 34:586–592
- 9) Castillo de Oyagüe R, Sánchez-Jorge MI, Sánchez Turrión A, Monticelli F, Toledano M, Osorio R (2009) Influence of CAM vs. CAD/CAM scanning methods and finish line of tooth preparation in the vertical misfit of zirconia bridge structures. *Am J Dent* 22:79–83
- 10) Sailer I, Feher A, Filser F, Gauckler L, Lüthy H, Hämmerle CH (2007) Five-year clinical results of zirconia frameworks for posterior fixed partial dentures. *Int J Prosthodont* 20:383–388
- 11) Lang NP, Kiel RA, Anderhalden K (1983) Clinical and microbiological effects of subgingival restorations with overhanging or clinically perfect margins. *J Clin Periodontol* 10:563–578
- 12) Koernschild KL, Campbell SD (2000) Periodontal tissue responses after insertion of artificial crowns and fixed partial dentures. *J Prosthet Dent* 84:492–498
- 13) Felton DA, Kanoy BE, Bayne SC, Wirthman GP (1991) Effect of in vivo crown margin discrepancies on periodontal health. *J Prosthet Dent* 65:357–364
- 14) Goldman M, Laosonthorn P, White RR (1992) Microleakage—full crowns and the dental pulp. *J Endod* 18:473–475
- 15) Rekow D, Thompson VP (2005) Near-surface damage—a persistent problem in crowns obtained by computer-aided design and manufacturing. *Proc Inst Mech Eng* 219:233–243
- 16) Molin MK, Karlsson SL (2008) Five-year clinical prospective evaluation of zirconia-based Denzir 3-unit FPDs. *Int J Prosthodont* 21:223–227
- 17) Tinschert J, Schulze KA, Natt G, Latzke P, Heussen N, Spiekermann H (2008) Clinical behavior of zirconia-based fixed partial dentures made of DC-Zirkon: 3-year results. *Int J Prosthodont* 21:217–222
- 18) Vult von Steyern P, Carlson P, Nilner K (2005) All-ceramic fixed partial dentures designed according to the DC-Zirkon technique. A 2-year clinical study. *J Oral Rehabil* 32:180–187
- 19) Reich S, Kappe K, Teschner H, Schmitt J (2008) Clinical fit of four-unit zirconia posterior fixed dental prostheses. *Eur J Oral Sci* 116:579–584

- 20) Beuer F, Aggstaller H, Edelhoff D, Gernet W, Sorensen J (2009) Marginal and internal fits of fixed dental prostheses zirconia retainers. *Dent Mater* 25:94–102
- 21) Beuer F, Naumann M, Gernet W, Sorensen JA (2008) Precision of fit: zirconia three-unit fixed dental prostheses. *Clin Oral Investig* 13:343–349
- 22) Tinschert J, Natt G, Mautsch W, Spiekermann H, Anusavice KJ (2001) Marginal fit of alumina- and zirconia-based fixed partial dentures produced by a CAD/CAM system. *Oper Dent* 26:367–374
- 23) Reich S, Wichmann M, Nkenke E, Proeschel P (2005) Clinical fit of all-ceramic three-unit fixed partial dentures, generated with three different CAD/CAM systems. *Eur J Oral Sci* 113:174–179
- 24) Wettstein F, Sailer I, Roos M, Hämmerle CHF (2008) Clinical study of the internal gaps of zirconia and metal frameworks for fixed partial dentures. *Eur J Oral Sci* 116:272–279
- 25) Bindl A, Mörmann WH (2007) Fit of all-ceramic posterior fixed partial denture frameworks in vitro. *Int J Periodontics Restorative Dent* 27:567–575
- 26) Gonzalo E, Suárez MJ, Serrano B, Lozano JF (2008) Marginal fit of zirconia posterior fixed partial dentures. *Int J Prosthodont* 21:398–399
- 27) Coli P, Karlsson S (2004) Precision of a CAD/CAM technique for the production of zirconium dioxide copings. *Int J Prosthodont* 17:577–580
- 28) Kohorst P, Brinkmann H, Li J, Borchers L, Stiesch M (2009) Marginal accuracy of four-unit zirconia fixed dental prostheses fabricated using different computer-aided design/computer-aided manufacturing systems. *Eur J Oral Sci* 117:319–325
- 29) Vigolo P, Fonzi F (2008) An in vitro evaluation of fit of zirconium-oxide-based ceramic four-unit fixed partial dentures, generated with three different CAD/CAM systems, before and after porcelain firing cycles and after glaze cycles. *J Prosthodont* 17:621–626
- 30) Att W, Komine F, Gerds T, Strub JR (2009) Marginal adaptation of three different zirconium dioxide three-unit fixed dental prostheses. *J Prosthet Dent* 101:239–247
- 31) Holmes JR, Bayne SC, Holland GA, Sulik WD (1989) Considerations in measurement of marginal fit. *J Prosthet Dent* 62:405–408
- 32) Kunii J, Hotta Y, Tamaki Y, Ozawa A, Kobayashi Y, Fujishima A, Miyazaki T, Fujiwara T (2007) Effect of sintering on the marginal and internal fit of CAD/CAM-fabricated zirconia frameworks. *Dent Mater J* 26:820–826
- 33) Persson ASK, Odén A, Andersson M, Sandborgh-Englund G (2009) Digitization of simulated clinical dental impressions: virtual three-dimensional analysis of exactness. *Dent Mater* 25:929–936
- 34) Witkowski S, Komine F, Gerds T (2006) Marginal accuracy of titanium copings fabricated by casting and CAD/CAM techniques. *J Prosthet Dent* 96:47–52
- 35) Komine F, Gerds T, Witkowski S, Strub JR (2005) Influence of framework configuration on the marginal adaption of zirconium dioxide ceramic anterior four-unit frameworks. *Acta Odontol Scand* 63:361–366
- 36) Dittmer MP, Borchers L, Stiesch M, Kohorst P (2009) Stresses and distortions within fixed dental prostheses due to the veneering process. *Acta Biomater* 5:3231–3239
- 37) Boening KW, Wolf BH, Schmidt AE, Kastner K, Walter MH (2000) Clinical fit of Procera AllCeram crowns. *J Prosthet Dent* 84:419–424
- 38) Coli P, Karlsson S (2004) Fit of a new pressure-sintered zirconium dioxide coping. *Int J Prosthodont* 17:59–64

- 39) Bornemann G, Lemelson S, Luthardt RG (2002) Innovative method for the analysis of the internal 3D fitting accuracy of Cerec-3 crowns. *Int J Comput Dent* 5:177-182
- 40) Luthardt RG, Bornemann G, Lemelson S, Walter MH, Hüls A (2004) An innovative method for evaluation of the 3-D internal fit of CAD/CAM crowns fabricated after direct optical versus indirect laser scan digitizing. *Int J Prosthodont* 17:680-685
- 41) Harris IR, Wickens JL (1994) A comparison of the fit of sparkeroled titanium copings and cast gold alloy copings. *Int J Prosthodont* 7:348-355
- 42) Laurent M, Scheer P, Dejou J, Laborde G (2008) Clinical evaluation of the marginal fit of cast crowns—validation of the silicone replica method. *J Oral Rehabil* 35:116-122
- 43) Rahme HY, Tehini GE, Adib SM, Ardo AS, Rifai KT (2008) In vitro evaluation of the “replica technique” in the measurement of the fit of Procera crowns. *J Contemp Dent Pract* 9:25-32
- 44) Groten M, Axmann D, Pröbster L, Weber H (2000) Determination of the minimum number of marginal gap measurements required for practical in vitro testing. *J Prosthet Dent* 83:40-49

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